

Chemistry 1210
Quantitative Determination of a Two-Component System

Date: _____

Name: _____

Station #: _____

OBJECTIVE: The objective of this experiment is to quantitatively determine the composition of a two-component system using spectrophotometric methods.

PROCEDURE: As in the Chemistry 1210 lab manual, page 36 plus additional procedure as noted below.

OBSERVATIONS:
Describe Nickel, Cobalt, and Mixture solutions here.

Two-Component System

PROCEDURE:

1. With a partner, determine the Absorbance vs. Wavelength of Ni and Cobalt from 350-550 nm. Since the Spec 20 must be zeroed each time, have the nickel and the cobalt samples ready to be measured for each wavelength.
 - a) Choose the desired wavelength (360 nm).
 - b) Zero the instrument .
 - c) Measure the absorbance of the cobalt known and the nickel knowns at that wavelength.
 - d) Change the wavelength (380 nm).
 - e) Re-zero the instrument .
 - f) Again measure the absorbance of nickel and cobalt at the new wavelength.
 - g) Repeat the procedure every 20 nm until a wavelength of 600 has been measured.
 - h) Find the wavelength of maximum absorbance (λ_{max}) more accurately by finding regions of high absorbance and re-scanning them in steps of 5 nm.

Note: To convert transmittance to absorbance use the following formula:

$$A=2\cdot\log(\%T)$$

3. **Now work on your own**, and use only one instrument for rest of the readings:

Measuring the knowns:

- a) Make sure the instrument is zeroed at the nickel λ_{\max} .
- b) At the nickel λ_{\max} , **determine % transmittance and/or the absorbance (see data section) of the known nickel and known cobalt twice each, the second of each with a fresh sample.**
- c) Re-zero the machine at the cobalt λ_{\max} .
- d) At the cobalt λ_{\max} , measure the **known** nickel and the **known** cobalt **as above**.

Measuring the unknowns:

- e) Empty the cuvettes containing known nickel and cobalt solutions and refill them (with correct rinsing) with your **unknown** cobalt. Also fill another cuvette with the unknown mixture.
- f) Since the instrument is still at the cobalt λ_{\max} , rezero , and measure the **unknown** cobalt and **unknown** mixture at the cobalt λ_{\max} . (Be sure to put your data in the correct table). **Take two readings of each as previously described.**
- g) At the nickel λ_{\max} , zero the instrument.
- h) Now measure the **unknown** mixture at that wavelength. Again, take two readings of each.
- i) **All cuvettes must be rinsed out thoroughly with distilled water and turned upside down in the test tube rack to indicate they are clean.**

DATA:

Determination of λ_{\max}

Use the extra space at the end to determine the two wavelengths more accurately by finding regions of highest absorbance for each metal and re-scanning them in 5 nm increments.

Wavelength	%T	Absorbance of known Co solution	%T	Absorbance of known Ni solution
360				
380				
400				
420				
440				
460				
480				
500				
520				
540				
560				
580				
600				

Attach a graph of absorbance vs. wavelength for both Co and Ni. Label each λ_{\max} clearly.

DATA:

Concentration of known nickel solution	
Concentration of known cobalt solution	

Ni λ_{\max}	
Co λ_{\max}	

*When less than 3 sig figs are obtained for absorbance, record %T, then calculate absorbance.
Note: Calculate average absorbance after converting %T to absorbance.*

% T or Absorbance of known nickel at the nickel λ_{\max}			Average absorbance: __
% T or Absorbance of known nickel at the cobalt λ_{\max}			Average absorbance: __

% T or Absorbance of known cobalt at the nickel λ_{\max}			Average absorbance: __
% T or Absorbance of known cobalt at the cobalt λ_{\max}			Average absorbance: __

Unknown# (Cobalt): _____

% T or Absorbance of unknown cobalt # _____ at the Cobalt λ_{\max}			Average absorbance: __
--	--	--	------------------------

Unknown# (mixture): _____

% T or Absorbance of unknown mixture at Nickel λ_{\max}			Average absorbance: __
% T or Absorbance of unknown mixture at Cobalt λ_{\max}			Average absorbance: __

CALCULATIONS:

3. In the space below, calculate the concentration of your unknown cobalt solution:

4. In the space below, calculate the concentration of the cobalt and of the nickel in your mixture.

RESULTS:

Unknown	Concentration	
Co# _____		
Mixture# _____	Ni	Co

DISCUSSION:

Give one source of error (beyond your reasonable control) in this experiment, and state how it would affect your results.